



Recombinant Human Neuralized-like protein 1A (NEURL)

Product Code	CSB-MP015721HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O76050
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	GNNFSSIPS LPRGNPSRAP RGHPQNLKDS IGGPFPVTSH RCHHKQKHCP AVLPSSGGLPA TPLLFHPHTK GSQILMDLSH KAVKRQASFC NAITFSNRPV LIYEQVRLKI TTKQCCWSGA LRLGFTSKDP SRIHPDSLKP YACPDLSVSQS GFWAKALPEE FANEGNIIAF WVDKKGRVFH RINDSAVMLF FSGVRTADPL WALVDVYGLT RGVQLLDSEL VLPDCLRPRS FTALRRPSLR READDARLSV SLCDLNVPGA DGDEAAPAAG CIPQNSLNS QHSRALPAQL DGDRLFHALR AGAHVRILDE QTVARVEHGR DERALVFTSR PVRVAETIFV KVTRSGGARP GALSFGVTTC DPGTLRPADL PFSPEALVDR KEFWAVCRVP GPLHSGDILG LVVNADGELH LSHNGAAAGM QLCVDASQPL WMLFGLHGTI TQIRILGSTI LAERGIPSLP CSPASTPTSP SALGSRLSDP LLSTCSSGPL GSSAGGTAPN SPVSLPESPV TPGLGQWSDE CTICYEHAVD TVIYTCGHMC LCYACGLRLK KALHACCPIC RRIKDIKT YRSS
Source	Mammalian cell
Target Names	NEURL1
Protein Names	Recommended name: Neuralized-like protein 1A Short name= h-neu Short name= h-neuralized 1 EC= 6.3.2.- Alternative name(s): RING finger protein 67
Expression Region	2-574
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.