



Recombinant Human E3 ubiquitin-protein ligase Praja-2 (PJA2)

Product Code	CSB-EP018050HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43164
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>SQYTEKEPA AMDQESGKAV WPKPAGGYQT ITGRRYGRRH AYVSFKPCMT RHERSLGRAG DDYEVLLEDD VPKENSSGSS PLDQVDSSLP SEPIFEKSET EIPTCGSALN QTTESSQS FV AVHHSEEGRD TLGSSTNLHN HSEGEYIPGA CSASSVQNGI ALVHTDSYDP DGKHGEDNDH LQLSAEVVEG SRYQESLGNT VFELENREAE AYTGLSPPVP SFNCEVRDEF EELDSVPLVK SSAGDTEFVH QNSQEIQRSS QDEMVSTKQQ NNTSQERQTE HSPEDAACGP GHICSEQNTN DREKNHGSSP EQVVRPKVRK LISSSQVDQE TGFNRHEAKQ RSVQRWREAL EVEESGSDDL LIKCEEYDGE HDCMFLDPY SRVITQRETE NNQMTSESGA TAGRQVDNT FWNGCGDYYQ LYDKDEDSSE CSDGEWSASL PHRFSGTEKD QSSSDESWET LPGKDENEPE LQSDSSGPEE ENQELSLQEG EQTSLEEGERI PWLQYNEVNE SSSDEGNEPA NEFAQPAFML DGNNNLEDDS SVSEDLVDW SLFDGFADGL GVAEAI SYVD PQFLTYMALE ERLAQAMETA LAHLESLAVD VEVANPPASK ESIDGLPETL VLEDHTAIGQ EQCCPICCSE YIKDDIATEL PCHHFFHKPC VSIWLQKSGT CPVCRRHFPP AVIEASAAPS SEPDPDAPPS NDSIAEAP</p>
Source	E.coli
Target Names	PJA2
Protein Names	Recommended name: E3 ubiquitin-protein ligase Praja-2 Short name= Praja2 EC= 6.3.2.- Alternative name(s): RING finger protein 131
Expression Region	2-708
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.