



Recombinant Escherichia coli Beta-barrel assembly-enhancing protease (bepA)

Product Code	CSB-MP355436ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P66948
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	DSA DTLPDMGTSA GSTLSIGQEM QMGDYYVRQL RGSAPLINDP LLTQYINSLG MRLVSHANSV KTPFHFFLIN NDEINAF AFF GGNVVLHSAL FRYSDNESQL ASVMAHEISH VTQRHLARAM EDQQRSAPLT WVGALGSILL AMASPQAGMA ALTGTLAGTR QGMISFTQQN EQEADRIGIQ VLQRSGFDPQ AMPTFLEKLL DQARYSSRPP EILLTHPLPE SRLADARNRA NQMRPMVVQS SEDFYLAKAR TLGMYNSGRN QLTSDLLDEW AKGNVRQQRA AQYGRALQAM EANKYDEARK TLQPLLAEP GNAWYLDLAT DIDLGQNKAN EAINRLKNAR DLRTNPVLQL NLANAYLQGG QPQEAANILN RYTFNNKDDS NGWDLAQAE AALNNRDQEL AARAEGYALA GRLDQAISLL SSASSQVKLG SLQQARYDAR IDQLRQLQER FKPYTKM
Source	Mammalian cell
Target Names	bepA
Protein Names	Recommended name: TPR repeat-containing protein yfgC
Expression Region	28-487
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.