



Recombinant Human Krueppel-like factor 4 (KLF4)

Product Code	CSB-BP012394HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43474
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MRQPPGESDM AVSDALLPSF STFASGPAGR EKTLRQAGAP NNRWREELSH MKRLPPVLPG RPYDLAAATV ATDLESGGAG AACGGSNLAP LPRRETEEFN DLLDLDFILS NSLTHPPESV AATVSSSSASA SSSSSPSSSSG PASAPSTCSF TYPIRAGNDP GVAPGGTGGG LLYGRESAPP PTAPFNLADI NDVSPSGGFV AELLRPELDP VYIPPQQPQP PGGGLMGKFV LKASLSAPGS EYGSPSVISV SKGSPDGSHV VVVAPYNGGP PRTCPIKQE AVSSCTHLGA GPPLSNGHRP AAHDFPLGRQ LPSRTTPTLG LEEVLSSRDC HPALPLPPGF HPHPGPNYPS FLPDQMQPQV PPLHYQGQSR GFVARAGEPC VCWPHFGTHG MMLTPSSPL ELMPPGSCMP EEPKPKRGRR SWPRKRTATH TCDYAGCGKT YTKSSHLKAH LRTHTGEKPY HCDWDGCGWK FARSEDELTRH YRKHTGHRPF QCQKCDRAFS RSDHLALHMK RHF</p>
Source	Baculovirus
Target Names	KLF4
Protein Names	Recommended name: Krueppel-like factor 4 Alternative name(s): Epithelial zinc finger protein EZF Gut-enriched krueppel-like factor
Expression Region	1-513
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.