



Recombinant Transforming protein rel polyprotein (V-REL)

Product Code	CSB-MP355608ASW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01126
Product Type	Recombinant Protein
Immunogen Species	Avian reticuloendotheliosis virus
Purity	>85% (SDS-PAGE)
Sequence	MDFLTNLRFTEG
Source	Mammalian cell
Target Names	V-REL
Protein Names	Recommended name: Transforming protein rel polyprotein Alternative name(s): p58 V-rel Cleaved into the following 3 chains: 1. Env polyprotein N-terminal 2. Transforming protein rel 3. Env polyprotein C-terminal
Expression Region	1-12
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.