



Recombinant Escherichia coli Transcriptional activator protein lysR (lysR)

Product Code	CSB-YP355928ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03030
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MAAVNLRHIE IFHAVMTAGS LTEAAHLLHT SQPTVSRELA RFEKVIKGLKLFERVRGRLHP TVQGLRLFEE VQRSWYGLDR IVSAAESLRE FRQGELSIACLPVFSQSFLP QLLQPFLARY PDVSLNIVPQ ESPLLEEWLS AQRHDLGLTE TLHTPAGTER TELLSLDEVCLPVPVGHPLAV KKVLTTPDDFQ GENYISLSRTDSYRQLLDQL FTEHQVKRRM IVETHSAASV CAMVRAGVGI SVVNPLTALDYAASGLVRR FSIAVPFTVS LIRPLHRPSS ALVQAFSGHL QAGLPKLVTS LDAILSSATT A
Source	Yeast
Target Names	lysR
Protein Names	Recommended name: Transcriptional activator protein lysR
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.