



Recombinant Hamster polyomavirus Major capsid protein VP1

Product Code	CSB-YP355947HBZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03092
Product Type	Recombinant Protein
Immunogen Species	Hamster polyomavirus (HaPyV) (Mesocricetus auratus polyomavirus 1)
Purity	>85% (SDS-PAGE)
Sequence	MCKPLWKPCP KPANVPKLIM RGGVGVLDLV TGEDSITQIE AYLNPRMGQN KPGTGTGQY YGFSQSIKVN SSSLTADEVKA NQLPYYSMAK IQLPTLNEDL TCDTLQMWEA VSVKTEVVG V GSELLNVHGYG SRSETKDIGI SKPVEGTTYH MFAVGGEPLD LQGLVQNYNA NYEAAIVSIK TVTGKAMTST NQVLDPTAKA KLDKDGRIPI EIWGPDP SKN ENSRYYG NFT GGTGTPPMQ FTNTLTTLVLL DENG VGPLCK GDGLYLSAAD VMGWYIEYNS AGWHWRGLPR YFNVTLRKRW VKNPYPVTSL LASLYNNMLP TIEGQPMERGE AAQVEEVRIY EGTEAVPGDP DVNRFIDKYG QQHTKPPAKP AN
Source	Yeast
Protein Names	Recommended name: Major capsid protein VP1 Alternative name(s): Major structural protein VP1
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.