



Recombinant Hepatitis B virus genotype D subtype adw Capsid protein (C)

Product Code	CSB-EP355964HEK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03147
Product Type	Recombinant Protein
Immunogen Species	Hepatitis B virus genotype D subtype adw (isolate United Kingdom/adyw/1979) (HBV-D)
Purity	>85% (SDS-PAGE)
Sequence	MDIDPYKEFG ATVELLSFLP SDFFPSVRDL LDAAAALYRD ALESPEHCSP HHTALRQAIL CWGDLMTLAT WVGTLNLEDP SRDLVVS YVN TNVGLKFRQL LWFHISCLTF GRETVLE YLV SFGVWIRTPP AYRPPNAPIL STLPETT VVR RRGRSPRRRT PSPRRRRS QS PRRRRRSQSRE SQC
Source	E.coli
Target Names	C
Protein Names	Recommended name: Capsid protein Alternative name(s): Core antigen Core protein HBcAg p21.5
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.