



Recombinant Friend spleen focus-forming virus Gag polyprotein (gag)

Product Code	CSB-BP356016FDD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03331
Product Type	Recombinant Protein
Immunogen Species	Friend spleen focus-forming virus (isolate 502) (FSFFV)
Purity	>85% (SDS-PAGE)
Sequence	GQTVTTPLS LTLEHWEDVQ RTASNQSVDV KKRRWVTFCS AEWPTFGVGVW PQDGTFNLDI ILQVKSQVFS PGPHGHPDQV PYIVTWEAIA YEPWPVKPF VSPKLSPTSPT APILPSGPST QPPRSALYP ALTPSIKPGP SPIMADLSLT FSQKTLRRTTE DRDRPPLTEM ATEKRPPPLL RFLPPLP
Source	Baculovirus
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Alternative name(s): Core polyprotein Cleaved into the following 2 chains: 1. Matrix protein p15 Short name= 2. MA 3. RNA-binding phosphoprotein p12 Alternative name(s): pp12
Expression Region	2-187
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.