



Recombinant Friend spleen focus-forming virus Gag polyprotein (gag)

| | |
|--------------------------|---|
| Product Code | CSB-EP356016FDD |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P03331 |
| Product Type | Recombinant Protein |
| Immunogen Species | Friend spleen focus-forming virus (isolate 502) (FSFFV) |
| Purity | >85% (SDS-PAGE) |
| Sequence | GQTVTTPLS LTLEHWEDVQ RTASNQSVDV KKRRWVTFCS AEWPTFGVGV PQDGTFNLDI ILQVKSQVFS PGPHGHPDQV PYIVTWEAIA YEPWPVKPF VSPKLSPTS APILPSGPST QPPRSALYP ALTPSIKPGP SPIMADLSLT FSQKTLRTE DRDRPPLTEM ATEKRPPPLL RFLPPLP |
| Source | E.coli |
| Target Names | gag |
| Protein Names | Recommended name: Gag polyprotein Alternative name(s): Core polyprotein Cleaved into the following 2 chains: 1. Matrix protein p15 Short name= 2. MA 3. RNA-binding phosphoprotein p12 Alternative name(s): pp12 |
| Expression Region | 2-187 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |