



Recombinant Moloney murine leukemia virus Gag polyprotein (gag), partial

Product Code	CSB-MP356017MHK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03332
Product Type	Recombinant Protein
Immunogen Species	Moloney murine leukemia virus (isolate Shinnick) (MoMLV)
Purity	≥85% (SDS-PAGE)
Sequence	PLRAGGNGQLQYWPFSSSDLYNWKNNNPSFSEDPGKLTALIESVLITHQPTW DCCQQLLG TLLTGEEKQRVLLLEARKAVRGDDGRPTQLPNEVDAAFPLERPDWDYTTQAGR NHLVHYRQ LLLAGLQNAGRSPTNLAKVKGITQGPNESPSAFLERLKEAYRRYTPYDPEDPG QETNVSM SFIWQSAPDIGRKLRLLEDLKNKTLGDLVREAEKIFNKRETPEEREERIRRETE EKEERR RTEDEQKEKERDRRRHREMSKLL
Source	Mammalian cell
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Short name= Pr65gagAlternative name(s): Core polyproteinCleaved into the following 4 chains: 1. Matrix protein p15 Short name= 2. MA 3. RNA-binding phosphoprotein p12Alternative n
Expression Region	216-478aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.