



# Recombinant Woolly monkey sarcoma virus Pol polyprotein (pol)

<b>Product Code</b>	CSB-YP356024WBA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03359
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Woolly monkey sarcoma virus (WMSV) (Smian sarcoma-associated virus)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	GSLGRGDRPG VYWEVDFTEV KPGRYGNRYL LVFIDTFSGW VEAFTPKTET ALTVCKRNST PLRIPKVLGS DNGPAFVAQV SQGLATQLGI NWKLHCAYRP QSSGQVERMN RTIKETLTKL ALETGGKDWV ALLPLALLRA KNTPSRFGLT PYEILYGGPP PILESGGTLG PDDNPLPVLV THLKALEVVR TQIWDQIKEV YKPGTVAIPH PFQVGDQVLV RRHRPGSLEP RWKGPYLVLV TTPTAVKVDG IAAWVHASHL KPAPPSAPDE SWELEKADHP LKLRIARRRN ESAK
<b>Source</b>	Yeast
<b>Target Names</b>	pol
<b>Protein Names</b>	Recommended name: Pol polyprotein Cleaved into the following 2 chains: 1. Reverse transcriptase/ribonuclease H Short name= 2. RT EC= 3. 2.7.7.49 EC= 4. 2.7.7.7 EC= 5. 3.1.26.4 6. Integrase Short name= 7. IN
<b>Expression Region</b>	1-294
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.