



Recombinant Bovine leukemia virus Putative uncharacterized protein PXBL-I

Product Code	CSB-BP356038BKI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03412
Product Type	Recombinant Protein
Immunogen Species	Bovine leukemia virus (isolate Japanese BLV-1) (BLV)
Purity	>85% (SDS-PAGE)
Sequence	ASVVGWGPS LHACPALVLS NDVTIDAWCP LCGPHERLQF ERIDTTLTCE THRINWTADG RPCGLNGTLF PRLHVSETRP QGPRRLWINC PLPAVRAQPG PVSLSPFERS PFQPYQCQLP SASSDGCPII GHGLLPWNNL VTHPVLGKVL ILNQMANFSL LPSFDTLLVD PLRLSVFAPD TRGAIRYLST LLTLCPATCI LPLGEPFSPN VPICRFPDRS NEPLSEFEL PPIQTPGLSW SVPAIDLFLT GPPSPCDRLH VWSSPQALQR FLHDPTLTWS ELVASRKIRL DSPLKLQLE NEWLSRLF
Source	Baculovirus
Protein Names	Recommended name: Putative uncharacterized protein PXBL-I
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.