



Recombinant Bovine leukemia virus Putative uncharacterized protein PXBL-III

Product Code	CSB-BP356039BKI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03414
Product Type	Recombinant Protein
Immunogen Species	Bovine leukemia virus (isolate Japanese BLV-1) (BLV)
Purity	≥85% (SDS-PAGE)
Sequence	HHSLCAFWVR IFPTQLPLFT PCCTPFLEIP ENLSSHPEEG CGSEVKIARA ATSLSFYSTL ARPRVLSPLT EVQNFLYKGM LGSKCAQYLF QKVL MNV FPC NKPQQRHSSH IQQHLGRLF
Source	Baculovirus
Protein Names	Recommended name: Putative uncharacterized protein PXBL-III
Expression Region	1-119
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.