



Recombinant Enterobacteria phage G4 Capsid protein (F)

Product Code	CSB-EP356107ECL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03642
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage G4 (Bacteriophage G4)
Purity	>85% (SDS-PAGE)
Sequence	SNVQTSADR VPHDLSHLVF EAGKIGRLKT ISWTPVVAGD SFECDMVGAI RLSPLRRGLA VDSRVDIFSF YIPHRHIYQG QWINFMKDGV NASPLPPVTC SSGWDSAAYL GTIPSSTLKV PKFLHQGYLN IYNNYFKPPW SDDLTYANPS NMPSEDYKWG VRVANLKSIW TAPLPPDTRT SENMTTGTST IDIMGLQAAY AKLHTEQERD YFMTRYRDIM KEFGGHTSYD GDNRPLLLMR SEFWASGYDV DGTDQSSLGQ FSGRVQQT FN HKVPRFYVPE HGVIMTLAVT RFPPTHEMEM HYLVGKENLT YTDIACDPAL MANLPPREVS LKEFFHSSPD SAKFKIAEGQ WYRTQPDRVA FPYNALDGFP FYSALPSTEL KDRVLVNTNN YDEIFQSMQL AHWNMQTKFN INVYRHMPPT RDSIMTS
Source	E.coli
Target Names	F
Protein Names	Recommended name: Capsid protein Alternative name(s): F protein GPF
Expression Region	2-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.