



# Recombinant Enterobacteria phage phiK Minor spike protein (H)

<b>Product Code</b>	CSB-EP356109EDM-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03649
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Enterobacteria phage phiK (Bacteriophage phi-K)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLGSIIGGIG SLLGGIASG GISSLLNKMFK SKMPEHAASS AGLTNGQGTI GMDTDAGIQS AIQGSNVPPA GQLPASNTSG VMADAGNMIR NAGRALLDGT IQAGSDKVKE ALIDKLVGGN DAKDRGKATR DYLAAPPEL NPWERAGAGA STAGLESSAQ NQQKEMLMQ LDNQKDIAM QMNNNLQIAG IQSATSQRNT KDSVYAQNEM LQYNQRESQA RVASILANTD LTTKQATHEI MRMALTRAQE TGQHLLNSQI MALEKKVYAE IGKIHQDTQN SRYGSSQVTA AAKDVTNMIT DASSGLADWA SQQWDSFFKD GKSDGIPLNT RK
<b>Source</b>	E.coli
<b>Target Names</b>	H
<b>Protein Names</b>	Recommended name: Minor spike protein Alternative name(s): H protein Pilot protein
<b>Expression Region</b>	1-332
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.