



Recombinant Human Sjogren syndrome/scleroderma autoantigen 1 (SSSCA1)

Product Code	CSB-BP022722HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60232
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ALNGAEVDD FSWEPPTEAE TKVLQARRER QDRISRLMGD YLLRGYRMLG ETCADCGTIL LQDKQRKIYC VACQELSDV DKDNPALNAQ AALSQAREHQ LASASELPLG SRPAPQPPVP RPEHCEGAAA GLKAAQGPPA PAVPPNTDVM ACTQTALLQK LTWASAEELGS STSLETSIQL CGLIRACAEA LRSLQQLQH
Source	Baculovirus
Target Names	SSSCA1
Protein Names	Recommended name: Sjogren syndrome/scleroderma autoantigen 1 Alternative name(s): Autoantigen p27
Expression Region	2-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren s syndrome. Subcellular localization has not yet been established.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.