



Recombinant Enterobacteria phage T4 Tail fiber protein Gp35 (35)

Product Code	CSB-YP356133EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03742
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MEKFMAEFGQ GYVQTPFLSE SNSVRYKISI AGSCPLSTAG PSYVKFQDNP VGSQTF SAGL HLRVFD PSTG ALVDSKSYAF STSNDTTSAA FVSFMNSLTN NRIVAILTSG KVNFPPEVVS WLRTAGTSAF PSDSILSRFD VSYAAFYTSS KRAIALEHVK LSNRKSTDDY QTILDVVFDS LEDVGATGFP RGTYESVEQF MSAVGGTND E IARLPTSAAI SKLSDYNLIP GDVLYLKAQL YADADLLALG TTNISIRFYN ASNGYISSTQ AEFTGQAGSW ELKEDYVVVP ENAVGFTIYA QRTAQAGQGG MRNLSFSEVS RGGISKPAE FGVNGIRVNY ICESASPPDI MVLPTQASSK TGKVFGEFR EV
Source	Yeast
Target Names	35
Protein Names	Recommended name: Tail fiber protein Gp35
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.