



Recombinant Enterobacteria phage T4 DNA replication protein (69)

Product Code	CSB-YP356279EDZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04523
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MAHFNECAHL IEGVDKAQNE YWDILGDEKD PLQVMLDMQR FLQIRLANVR EYCYHPDKLE TAGDVVSWMR EQKDCIDDEF RELTSLGEM SRGEKEASAV WKKWKARYIE AQEKRIDEMS PEDQLEIKFE LVDIFHFVLN MFVGLGMNAE EIFKLYYLKN KHILNVKIMD IKQKFYRTYI VKVRTPKG VF WYAGKHESFI VNPYNDKYPG SGKILWNIYR KYGFNYKIRW SKCHGSREKS YEVERELISA LKRKHPDTCI NISPGGQGGG GRKWTEQQRL EHKLRLLNPE TKTRMKNSQR IAQNRAERKA RQSEVMKKFY SNGGNKKISE GTSRAQRKAP HWHEPLKSEI HELWVSLGKP ATGPVVKALK GKYDVTSSAL KNLIYLF RKE DV
Source	Yeast
Target Names	69
Protein Names	Recommended name: DNA replication protein Alternative name(s): Protein Gp69
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.