



Recombinant Human PRAME family member 2 (PRAMEF2)

Product Code	CSB-MP018615HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60811
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSIQAPPRLL ELAGQSLLRD QALSISAMEE LPRVLYLPLF REAFSRRHFQ TLTVMVQAWP FTCLPLVSLM KTLHLEPLKA LLEGLHMLLT QKDRPRRWKL QVLDLRDVDE NFWARWPGAW ALSCFPEAMS KRQTAEDCPR TGEHQPLKVF IDICLKEIPQ DECLRYLFQW VYQRRGLVHL CCSKLVNYLT PIKYLRKSLK IYINSIGEL EIHNTCWP HL IRKLYCYLKE MKTLCKLVFS RCHHYTSDNE LEGWLVTRFT SVFLRLEHLQ LLKIKLITFF SGHLEQLIRC LQNPLENLEL TCGNLEEDL KCLSQFPSLG YLKHNLNSYV LLFRISLEPL GALLEKIAAS LETLVLEG CQ IHYSQLSAIL PGLSCCSQLT TFYFGSNCMS IDALKDLLRH TSGLSKLSLE TYPAPEESLN SLVRVNWEIF TPLRAELMCT LREFRQPKRI FIGPTPCPSC GSSPSELEL HLCC
Source	Mammalian cell
Target Names	PRAMEF2
Protein Names	Recommended name: PRAME family member 2
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.