




Recombinant Influenza A virus Hemagglutinin (HA)

Product Code	CSB-YP356317IDG
Relevance	Binds to sialic acid-containing receptors on the cell surface, bringing about the attachment of the virus particle to the cell. This attachment induces virion internalization of about two third of the virus particles through clathrin-dependent endocytosis and about one third through a clathrin- and caveolin-independent pathway. Plays a major role in the determination of host range restriction and virulence. Class I viral fusion protein. Responsible for penetration of the virus into the cell cytoplasm by mediating the fusion of the membrane of the endocytosed virus particle with the endosomal membrane. Low pH in endosomes induces an irreversible conformational change in HA2, releasing the fusion hydrophobic peptide. Several trimers are required to form a competent fusion pore.
Abbreviation	Recombinant Influenza A virus Hemagglutinin protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04664
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/England/878/1969 H3N2)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QDLPGNDNSTATLCLGHHAVPNGTLVKTITNDQIEVTNATELVQSSSTGKICNN PHRILDGINCTLIDALLGDPHCDVFQDETWDLFVERSKAFSNCYPYDVPDYASL RSLVASSGTLEFITEGFTWTGVTQNGGSNACKRGPDSGFFSRLNWLTKSGST YPVLNVTMPNNDNFDKLYIWGVHHPSTNQEQTSLYVQASGRVTVSTRRSQQT IIPNIGSRPWVRLSSRISYWTIVKPGDVLVINSNGNLIAPRGYFKMRTGKSSIM RSDAPIDTCISECITPNGSIPNDKPFQNVNKITYGACPKYVKQNTLKLATGMRN VPEKQT
Research Area	Microbiology
Source	Yeast
Target Names	HA
Protein Names	Recommended name: Hemagglutinin Cleaved into the following chain: 1. Hemagglutinin HA1 chain
Expression Region	1-328aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	38.1 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.