



Recombinant Escherichia coli Lipoprotein 28 (nlpA)

Product Code	CSB-MP356377ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04846
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	CDQSSSD AKHIKVGVIN GAEQDVAEVA KKVAKKEYGL DVELVGFSGS LLPNDATNHG ELDANVFQHR PFLEQDNQAH GYKLVAVGNT FVFPMAGYSK KIKTVAQIKE GATVAIPNDP TNLGRALLL QKEKLITLKE GKGLLPTALD ITDNPRHLQI MELEGAQLPR VLDDPKVDVA IISTTYIQQT GLSPVHDSVF IEDKNSPYVN ILVAREDNKN AENVKEFLQS YQSPEVAKAA ETIFNGGAVP GW
Source	Mammalian cell
Target Names	nlpA
Protein Names	Recommended name: Lipoprotein 28
Expression Region	24-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.