



Recombinant Enterobacteria phage P22 Essential recombination function protein (erf)

Product Code	CSB-MP356398EDG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04892
Product Type	Recombinant Protein
Immunogen Species	Salmonella phage P22 (Bacteriophage P22)
Purity	≥85% (SDS-PAGE)
Sequence	MSKEFYARLA EIQEHLNAPK NQYNSFGKYK YRSCEDILEG VKPLLKGLFL SISDEIVLIG DRYVVKATAT ITDGENSHSA SAIAREEENK KGMDAAQVTG ATSSYARKYC LNGLFGIDDA KDADTEEHKQ QQNAAPAKQT KSSPSSPAPE QVLKAFSEYA ATETDKKLI ERYQHDWQLL TGHDDQTKC VQVMNIRINE LKQVA
Source	Mammalian cell
Target Names	erf
Protein Names	Recommended name: Essential recombination function protein
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.