



Recombinant Escherichia coli Exodeoxyribonuclease V alpha chain (recD)

Product Code	CSB-BP356429ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04993
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKLQKQLLEA VEHKQLRPLD VQFALTVAGD EHPAVTLAAA LLSHDAGEGH VCLPLSRLEN NEASHPLLAT CVSEIGELQN WEECLLASQA VSRGDEPTPM ILCGDRLYLN RMWCNERTVA RFFNEVNHAI EVDEALLAQT LDKLFPVSDE INWQKVAAAV ALTRRISVIS GPGGTGKTTT VAKLLAALIQ MADGERCRIR LAAPTGKAAA RLTESLGKAL RQLPLTDEQK KRIPEDASTL HRLLGAQPGS QRLRHHAGNP LHLVDLVVDE ASMIDLPMMMS RLIDALPDHA RVIFLGDRDQ LASVEAGAVL GDICAYANAG FTAERARQLS RLTGTHVPAG TGTEAASLRD SLCLLQKSYR FGSDSGIGQL AAAINRGDKT AVKTVFQQDF TDIEKRLLS GEDYIAMLEE ALAGYGRYLD LLQARAEPDL IIQAFNEYQL LCALREGPFG VAGLNERIEQ FMQQKRKIHR HPHSRWYEGR PVMIARNSA LGLFNGDIGI ALDRGQGTRV WFAMPDGNIK SVQPSRLPEH ETTWAMTVHK SQGSEFDHAA LILPSQRTPV VTRELVYTAV TRARRRLSLY ADERILSAAI ATRTERRSGL AALFSSRE
Source	Baculovirus
Target Names	recD
Protein Names	Recommended name: Exodeoxyribonuclease V alpha chain EC= 3.1.11.5 Alternative name(s): Exodeoxyribonuclease V 67 kDa polypeptide
Expression Region	1-608
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.