



Recombinant Escherichia coli Exodeoxyribonuclease I (sbcB)

Product Code	CSB-BP356430ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04995
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MMNDGKQQST FLFDYETFG THPALDRPAQ FAAIRTDSEF NVIGEPEVFY CKPADDYLPQ PGAVLITGIT PQEARAKGEN EAAFAARIHS LFTVPKTCIL GYNNVRFDDDE VTRNIFYRNF YDPYAWSWQH DNSRWDLDDV MRACYALRPE GINWPENDDG LPSFRLEHLT KANGIEHSNA HDAMADVYAT IAMAKLVKTR QPRLFDYLFT HRNKHKLML IDVPQMKPLV HVSGMFGAWR GNTSWVAPLA WHPENRNAVI MVDLAGDISP LLELSDTLR ERLYTAKTDL GDNAAVPVKL VHINKCPVLA QANTLRPEDA DRLGINRQHC LDNLKILREN PQVREKVVAI FAEAEPFTPS DNVDALYNG FFSDADRAAM KIVLETEPRN LPALDITFVD KRIEKLLFNY RARNFPGTLD YAEQQRWLEH RRQVFTPEFL QGYADELQML VQQYADDKEK VALLKALWQY AEEIV
Source	Baculovirus
Target Names	sbcB
Protein Names	Recommended name: Exodeoxyribonuclease I Short name= Exonuclease I EC= 3.1.11.1 Alternative name(s): DNA deoxyribophosphodiesterase Short name= dRPase
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.