



# Recombinant Transposon Tn7 transposition protein tnsC (tnsC)

<b>Product Code</b>	CSB-EP356642ENL-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P05846
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSATRIQAVY RDTGVEAYRD NPFIEALPPL QESVNSAASL KSSLQLTSSD LQKSRVIRAH TICRIPDDYF QPLGTHLLLS ERISVMIRGG YVGRNPKTGD LQKHLQNGYE RVQTGELETF RFEEARSTAQ SLLIGCSGS GKTTSLHRIL ATYPQVIYHR ELNVEQVVYL KIDCSHNGSL KEICLNFFRA LDRALGSNYE RRYGLKRHGI ETMLALMSQI ANAHALGLLV IDEIQHLSRS RSGGSQEMLN FFVTMVNIIG VPVMLIGTPK AREIFEADLR SARRGAGFGA IFWDPIQQTQ RGKPNQEWIA FTDNLWQLQL LQRKDALLSD EVRDVWYELS QGVMDIVVKL FVLAQLRALA LGNERITAGL LRQVYQDELK PVHPMLEALR SGIPERIARY SDLVVPEIDK RLIQLQLDIA AIQEQTPEEK ALQELDTEDQ RHLYLMLKED YDSSLLIPTI KKAFSQNPTM TRQKLLPLVL QWLMEGETVV SELEKPSKSK KVSIAIKVVKP SDWDSLPTD LRYIYSQRQP EKTMHRLKG KGVIVDMASL FKQAG
<b>Source</b>	E.coli
<b>Target Names</b>	tnsC
<b>Protein Names</b>	Recommended name: Transposon Tn7 transposition protein tnsC Short name= Protein E
<b>Expression Region</b>	1-555
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.