



Recombinant Human Microphthalmia-associated transcription factor (MITF)

Product Code	CSB-YP014595HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75030
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MLEMLEYNHYQVQTHLENPTKYHIQQAQRQQVKQYLSTTLANKHANQVLSLP CPNQPGDH VMPPVPGSSAPNSPMAMLTLSNCEKEGFYKFEEQNRAESECPGMNTHSRA SCMQMDDVI DDIISLESSYNEEILGLMDPALQMANTLPVSGNLIDLYGNQGLPPPGLTISNSCP ANLPN IKRELTESEARALAKERQKKDNHNLIERRRRFNINDRIKELGTLPKSNPDPMR WNKGTI LKASVDYIRKLQREQQRAKELENRQKKLEHANRHLLLRIQELEMQARAHGLSLI PSTGLC SPDLVNRIIKQEPVLENCSDLLQHHADLTCTTTLDLTDGTITFNNNLGTGTEA NQAYSV PTKMGSKLEDILMDDTLSPVGVTDPLLSSVSPGASKTSSRRSSMSMEETEHTC</p>
Source	Yeast
Target Names	MITF
Protein Names	Recommended name: Microphthalmia-associated transcription factor Alternative name(s): Class E basic helix-loop-helix protein 32 Short name= bHLHe32
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein of Isoform M2
Target Details	This gene encodes a transcription factor that contains both basic helix-loop-helix and leucine zipper structural features. It regulates the differentiation and development of melanocytes retinal pigment epithelium and is also responsible for pigment cell-specific transcription of the melanogenesis enzyme genes. Heterozygous mutations in the this gene cause auditory-pigmentary syndromes, such as Waardenburg syndrome type 2 and Tietz syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.