



Recombinant Rabies virus Nucleoprotein (N)

Product Code	CSB-BP356708RAG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06025
Product Type	Recombinant Protein
Immunogen Species	Rabies virus (strain Pasteur vaccins / PV) (RABV)
Purity	>85% (SDS-PAGE)
Sequence	MDADKIVFKV NNQVVSLKPE IVDQYEYKY PAIKDLKKPC ITLGKAPDLN KAYKSVLSCM SAAKLDPPDV CSYLAAAMQF FEGTCPEDWT SYGIVIARKG DKITPGSLVE IKRTDVEGNW ALTGGMELTR DPTVPEHASL VGLLLSLYRL SKISGQSTGN YKTNIADRIE QIFETAPFVK IVEHHTLMTT HKMCANWSTI PNFRFLAGTY DMFFSRIEHL YSAIRVGTVV TAYEDCSGLV SFTGFIKQIN LTAREAILYF FHKNFEEEEIR RMFEPGQETA VPHSYFIHFR SLGLSGKSPY SSNAVGHVFN LIHFVGCYMG QVRSLNATVI AACAPHEMSV LGGYLGEEFF GKGTFFERRFF RDEKELQEYE AAELTKTDVA LADDGTVNSD DEDYFSGETR SPEAVYTRII MNGGRLKRSH IRRYVSVSSN HQARPNSFAE FLNKTYSSDS
Source	Baculovirus
Target Names	N
Protein Names	Recommended name: Nucleoprotein Short name= NP Alternative name(s): Nucleocapsid protein Short name= Protein N
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.