



Recombinant *Saccharomyces cerevisiae* cAMP-dependent protein kinase type 2 (TPK2)

Product Code	CSB-MP356756SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06245
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEFVAERAQP VGQTIQQQNV NTYGQGVLPQ HHDLQQRQQQ QQQRQHQQLL TSQLPQKSLV SKGKYTLHDF QIMRTLGTGS FGRVHLVRSV HNGRYAIV LKKQQVVKMK QVEHTNDERR MLKLVEHPFL IRMWGTFQDA RNIFMVMYI EGGELFSLR KSQRFPNPVA KFYAAEVILA LEYLHAHNII YRDLKPENIL LDRNGHIKIT DFGFAKEVQT VTWTLCGTPD YIAPEVITTK PYNKSVDWWS LGVLIYEMLA GYTPFYDTP MKTYEKILQG KVVYPPYFHP DVVDLLSKLI TADLRRIGN LQSGSRDIKA HPWFSEVVWE RLLAKDIETP YEPPITSGIG DTSFLDQYPE EQLDYGIQGD DPYAEYFQDF
Source	Mammalian cell
Target Names	TPK2
Protein Names	Recommended name: cAMP-dependent protein kinase type 2 Short name= PKA 2 EC= 2.7.11.11
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.