



Recombinant Hordeum vulgare Serpin-Z4 (PAZ1)

Product Code	CSB-EP356767HWQ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06293
Product Type	Recombinant Protein
Immunogen Species	Hordeum vulgare (Barley)
Purity	>85% (SDS-PAGE)
Sequence	MATTLATDVR LSIHQTRFA LRLRSAISSN PERAAGNVAF SPLSLHVALS LITAGAAATR DQLVAILGDG GAGDAKELNA LAEQVVQFVL ANESSTGGPR IAFANGIFVD ASLSLKPSFE ELAVCQYKAK TQSVDFQHKT LEAVGQVNSW VEQVTTGLIK QILPPGSVDN TTKLILGNAL YFKGAWDQKF DESNTKCDSF HLLDGSSIQT QFMSSTKKQY ISSSDNLKVL KLPYAKGHDK RQFSMYILLP GAQDGLWSLA KRLSTEPEFI ENHIPKQTV E VGRFQLPKFK ISYQFEASSL LRALGLQLPF SEEADLSEMV DSSQGLEISH VFHKSFVEVN EEGTEAGAAT VAMGVAMSMP LKVDLVDFVA NHPFLFLIRE DIAGVVV FVG HVTNPLISA
Source	E.coli
Target Names	PAZ1
Protein Names	Recommended name: Serpin-Z4 Alternative name(s): BSZ4 HorvuZ4 Major endosperm albumin Protein Z4 Short name= Protein Z
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.