



# Recombinant Human CCR4-NOT transcription complex subunit 3 (CNOT3)

<b>Product Code</b>	CSB-EP005666HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	O75175
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MADKRKLQGE IDRCLKKVSE GVEQFEDIWQ KLNHAANANQ KEKYEADLKK EIKKLQRLRD QIKTWVASNE IKDKRQLIDN RKLIENTQMER FKVVERETKT KAYSKEGLGL AQKVDPAQKE KEEVGQWLTN TIDTLNMQVD QFESEVESLS VQTRKKKGDK DKQDRIEGLK RHIEKHRYHV RMLETILRML DNDSILVDAI RKIKDDVEYY VDSSQDPDFE ENEFLYDDLD LEDIPQALVA TSPPSHSHME DEIFNQSSST PTSTTSSSPI PPSPANCTTE NSEDDKKRGR STDSEVSQSP AKNGSKPVHS NQHPQSPAVP PTYPSGPPPA ASALSTTPGN NGVPAPAAPP SALGPKASPA PSHNSGTPAP YAQAVAPPAP SGPSTTQPRP PSVQPSGGGG GGSGGGGSSS SSNSAGGGA GKQNGATSYS SVVADSPAEV ALSSSGGNNNA SSQALGPPSG PHNPPPSTSK EPSAAAPTGA GGVAPGSGNN SGGPSLLVPL PVNPPSSPTP SFSDAKAAGA LLNGPPQFST APEIKAPEPL SSLKSMARA AISSGIEDPV PTLHLTERDI ILSSTSAPPA SAQPPLQLSE VNIPLSLGVC PLGPVPLTKE QLYQQAMEEA AWHHMPHPSD SERIRQYLPR NPCPTPPYHH QMPPPHSDTV EfyQRLSTET LFFIFYYLEG TKAQYLAACA LKKQSWRFHT KYMMWFQRHE EPKTITDEFE QGTYYFDYE KWGQRKKEGF TFEYRYLEDR DLQ
<b>Source</b>	E.coli
<b>Target Names</b>	CNOT3
<b>Protein Names</b>	Recommended name: CCR4-NOT transcription complex subunit 3 Alternative name(s): CCR4-associated factor 3 Leukocyte receptor cluster member 2
<b>Expression Region</b>	1-753
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.