



Recombinant Gag polyprotein (gag)

Product Code	CSB-EP356824ASS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06444
Product Type	Recombinant Protein
Immunogen Species	Avian myelocytomatosis virus HBI
Purity	>85% (SDS-PAGE)
Sequence	MEAVIKVISS ACKTYCGKTS PSKKEIGAML SLLQKEGLLM SPSDLYSPGS WDPITAALSQ RAMVLGKSGE LKTWGLVLGA LKAAREEQVT SEQAKFWLGL GGGRVSPPGP ESIEKPATER RIDKGEEVGE TTVQRDAKMA PEETATPKTV GTSCYHCGTA IGCNCATASA PPPPY
Source	E.coli
Target Names	gag
Protein Names	Recommended name: Gag polyprotein Cleaved into the following 3 chains: 1. Matrix protein p19 2. Core protein p10 3. Core protein p27
Expression Region	1-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.