



Recombinant Chlamydia trachomatis serovar L2 Major outer membrane porin (ompA)

Product Code	CSB-EP356881DSG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P06597
Product Type	Recombinant Protein
Immunogen Species	Chlamydia trachomatis serovar L2 (strain 434/Bu / ATCC VR-902B)
Purity	>85% (SDS-PAGE)
Sequence	LPVGNPAE PSLMIDGILW EGFGGDPCDP CTTWCDAISM RMGYYGDFVF DRVLQTDVNK EFQMGAKPTT ATGNAAAPST CTARENPAYG RHMQDAEMFT NAAYMALNIW DRFDVFCTLG ATSGYLKGNs ASFNLVGLFG DNENHATVSD SKLVPNMSLD QSVVELYTDt TFAWSAGARA ALWECGCATL GASFQYAQSK PKVEELNVLC NAAEFTINKP KGYVGQEFPL DLKAGTDGVT GTKDASIDYH EWQASLALSY RLMNFPTYIG VKWSRASFDA DTIRIAQPKS ATTVFDVTTL NPTIAGAGDV KASAEGQLGD TMQIVSLQLN KMKSRKSCGI AVGTTIVDAD KYAVTVETRL IDERAAHVNA QFRF
Source	E.coli
Target Names	ompA
Protein Names	Recommended name: Major outer membrane porin Short name= MOMP
Expression Region	23-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.