



Recombinant Human Sperm-associated antigen 7 (SPAG7)

Product Code	CSB-MP022468HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75391
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ADLLGSILS SMEKPPSLGD QETRRKAREQ AARLKKLQEQ EKQQKVEFRK RMEKEVSDFI QDSGQIKKKF QPMNKIERSI LHDVVEVAGL TSFSFGEDDD CRYVMIFKKE FAPSDEELDS YRRGEEWDPQ KAEKRLKE LAQRQEEAAA QQGPVVSPA SDYKDKYSHL IGKGAAKDAA HMLQANKTYG CVPVANKRDT RSIEEAMNEI RAKKRLRQSG EELPPTS
Source	Mammalian cell
Target Names	SPAG7
Protein Names	Recommended name: Sperm-associated antigen 7
Expression Region	2-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.