



Recombinant Enterobacteria phage D108 Repressor protein CI (CI)

Product Code	CSB-MP357004ECG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07040
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage D108 (Bacteriophage D108)
Purity	≥85% (SDS-PAGE)
Sequence	MMSFEIKEWF NAKELEGMPG VPKLATNITR KAVAEDWVKR QRHGGKGVAY EYHINSLPEE TRRAIKGASL SDKPVHTSIV HTVDERLIYA MSFLTPDEQA AAVEIIRVAG IKGLMPTIVS KDKALEALGI TVEQQKTLQT LQALPPEKVR EILSQYEGKE HNFVRENDV KKAIV
Source	Mammalian cell
Target Names	repc
Protein Names	Recommended name: Repressor protein CI
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.