



Recombinant Human Serpin B7 (SERPINB7)

Product Code	CSB-MP021075HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75635
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASLAAANAE FCFNLFREMD DNQGNGNVFF SSLSLFAALA LVRLGAQDDSL LSQIDKLLHV NTASGYGNSS NSQSGLQSQL KRVFSDINAS HKDYDLSIVN GLFAEKVYGF HKDYIECAEK LYDAKVERVD FTNHLEDTRR NINKWVENET HGKIKNVIGE GGISSSAVMV LVNAVYFKGK WQSAFTKSET INCHFKSPKC SGKAVAMMHQ ERKFNLVIE DPSMKILELR YNGGINMYVL LPENDLSEIE NKLTFQNLME WTNPRRMTSK YVEVFFPQFK IEKNYEMKQY LRALGLKDIF DESKADLSGI ASGGRLYISR MMHKSYIEVT EEGTEATAAT GSNIVEKQLP QSTLFRADHP FLFVIRKDDI ILFSGKVSCP
Source	Mammalian cell
Target Names	SERPINB7
Protein Names	Recommended name: Serpin B7 Alternative name(s): Megsin TP55
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.