



Recombinant Maltodextrin phosphorylase (malP)

Product Code	CSB-MP357022KBG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07094
Product Type	Recombinant Protein
Immunogen Species	Klebsiella pneumoniae
Purity	≥85% (SDS-PAGE)
Sequence	MSQTSFNKAQ FQAALTRQWQ HFGLQSASEM TQRQWWRAVS GALSELSAQ PVAKATEGER HVNYISMEFL IGRLTGNNLL NLGWYQQVSD ELQAHDVNLT DLLEEEID
Source	Mammalian cell
Target Names	malP
Protein Names	Recommended name: Maltodextrin phosphorylase EC= 2.4.1.1
Expression Region	1-108
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.