



Recombinant Human Keratin, type I cuticular Ha6 (KRT36)

| | |
|--------------------------|--|
| Product Code | CSB-EP012552HU-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | O76013 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATQTCTPTF STGSIKGLCG TAGGISRVSS IRSVGSCRVP SLAGAAGYIS SARSGLSGLG SCLPGSYLSS ECHTSGFVGS GWFCEGSFN GSEKETMQFL NDRLANYLEK VRQLERENAE LESRIQEWYE FQIPYICPDY QSYFKTIEDF QQKILLTKSE NARLVLQIDN AKLAADDFRT KYETELSLRQ LVEADINGLR RILDELTLCCK ADLEAQVESL KEELMCLKKN HEEEVSVLRC QLGDRLNVEV DAAPPVDLNK ILEDMRCQYE ALVENNRRDV EAWFNTQTEE LNQQVVSSSE QLQCCQTEII ELRRTVNALE IELQAQHSMR NSLESTLAET EARYSSQLAQ MQCLISNVEA QLSEIRCDLE RQNQEYQVLL DVKARLEGEI ATYRHLLIGE DCKLPPQPCA TACKPVIRVP SVPPVPCVPS VPCTPAPQVG TQIRTITEEI RDGKVISSRE HVQSRPL |
| Source | E.coli |
| Target Names | KRT36 |
| Protein Names | Recommended name: Keratin, type I cuticular Ha6 Alternative name(s): Hair keratin, type I Ha6 Keratin-36 Short name= K36 |
| Expression Region | 1-467 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Target Details | This protein is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.