



Recombinant Enterobacteria phage T2 Receptor-recognizing protein (38)

Product Code	CSB-BP357239EDX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07875
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T2 (Bacteriophage T2)
Purity	≥85% (SDS-PAGE)
Sequence	MAIVGVPGWI GESAVNETGQ RWMDAAMRAV HVSVPGWMSS MAGQSKEIYL SIGANHNYDR NSLINWMRAQ GGAPVVITIT GDLVSNSTGN ACLEFPSDLP NAYIQLIINS GVTYVGRGGN GSTNSSAGGN GGTAIHNAAG TKLRIRNNGA IAGGGGGGGA ASLKNSYPTN GSCGGGGGGRP FGVGGKIGSD SILSGSNASL TDAGTGGTTF QYGAGNGGNV GAGGGRGWGK NVYTSEGGAA GAAVTGNAPN WQNVGTIYGS RV
Source	Baculovirus
Target Names	38
Protein Names	Recommended name: Receptor-recognizing protein Alternative name(s): Protein Gp38
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.