



Recombinant Glucose 1-dehydrogenase B (gdhB)

Product Code	CSB-MP357272BQV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P07999
Product Type	Recombinant Protein
Immunogen Species	Bacillus megaterium
Purity	>85% (SDS-PAGE)
Sequence	MYKDLEGKVV VITGSSTGLG KSMAIRFATE KAKVVVNYRS KEDEANSVLE EEIKKVGGEA IAVKGDVTVE SDVINLVQSA IKEFGKLDVM INNAGMENPV SSHEMSLSDW NKVIDTNLTG AFLGSREAIK YFVENDIKGT VINMSSVHEW KIPWPLFVHY AASKGGMKLM TETLALEYAP KGIRVNNIGP GAINTPINAE KFADPEQRAD VESMIPMGYI GEPEEIAAVA WLASSEASYV TGITLFADGG MTQYPSFQAG RG
Source	Mammalian cell
Target Names	gdhB
Protein Names	Recommended name: Glucose 1-dehydrogenase B EC= 1.1.1.47
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.