



Recombinant Mobilization protein A (mobA)

Product Code	CSB-EP357300ENL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08098
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MGLSHKERRE YGRSYKRDNN RTAPFSCWPR WIRPGGWSGH GGMSGTTRTG VTARLAFWGG WLSPCLSRAG EKPSCVRTST RTQIITCCSR LTCSGLPGAR AGWLSGLCWR GSGVLLLDGH FAAADEAVNR KLNRRAGAQP ATGEMTDVRQ PAFRGAGTVN ALARFFHGRG PETAGVFLGK DEQGEPVLVP RDTWRKTNIQ ILGLPGSGKS VMATNALIRS VRDFGDAVVY FDQRGPVAAL LPAHCPNFTL LDLRPGKPAQ LNLFRDLQY ALKNLLVAGF NLSETGHVAD HYRISEQKAA KLIAEQVSAG GKHSAGAGGG VRAAGGPEKG REGADHQAGE RGRPERPADG QRHRRGGDNQ RRRLLYVIGS MDDEAVIRVQ KMLFARCAQI IARDEFRRW PHASIMLDEI KYLLSKYVLN ALGTVRSRDC NLRLAHHAGR LRAAGQDLPA DFKTTVLDN TPIRWFYRAA SQESRSGAGQ TGEIRVDVER RGPAGRRGTW SISAGTASSR RSRPLFDVNT LQHLLTGSR
Source	E.coli
Target Names	mobA
Protein Names	Recommended name: Mobilization protein A Short name= Protein B
Expression Region	1-529
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.