



Recombinant Dermatophagoides pteronyssinus Peptidase 1 (DERP1)

Product Code	CSB-YP357315DJH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08176
Product Type	Recombinant Protein
Immunogen Species	Dermatophagoides pteronyssinus (European house dust mite)
Purity	>85% (SDS-PAGE)
Sequence	TN ACSINGNAPA EIDLRQMRTV TPIRMQGGCG SCWAFSGVAA TESAYLAYRN QSLDLAEQEL VDCASQHGCH GDTIPRGIEY IQHNGVVQES YYRYVAREQS CRRPNAQRFG ISNYCQIYPP NVNKIREALA QTHSAIAVII GIKDLDAFRH YDGRITIIQRD NGYQPNYHAV NIVGYSNAQG VDYWIVRNSW DTNWGDNGYG YFAANIDLMM IEEYPYVVIL
Source	Yeast
Target Names	DERP1
Protein Names	Recommended name: Peptidase 1 EC= 3.4.22.65 Alternative name(s): Allergen Der p I Major mite fecal allergen Der p 1 Allergen= Der p 1
Expression Region	99-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.