



Recombinant *Pseudomonas aeruginosa* Ornithine carbamoyltransferase, catabolic (arcB)

Product Code	CSB-EP357346EZC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08308
Product Type	Recombinant Protein
Immunogen Species	<i>Pseudomonas aeruginosa</i> (strain ATCC 15692 / DSM 22644 / CIP 104116 / JCM 14847 / LMG 12228 / 1C / PRS 101 / PAO1)
Purity	>85% (SDS-PAGE)
Sequence	AFNMHNRNL LSLMHHSTRE LRYLLDLSRD LKRAKYTGTE QQHLKRKNIA LIFEKTSTRT RCAFEEVAAYD QGANVTYIDP NSSQIGHKES MKDTARVLGR MYDAIEYRGF KQEIVEELAK FAGVPVFNGL TDEYHPTQML ADVLTMREHS DKPLHDISYA YLGDARNNMG NSLLLIGAKL GMDVRIAAPK ALWPHDEFVA QCKKFAEESG AKLTLTEDPK EAVKGVDFVH TDVWVSMGEP VEAWGERIKE LLPYQVNMEI MKATGNPRAK FMHCLPAFHN SETKVGKQIA EQYPNLANGI EVTEDVFESP YNIAFEQAEN RMHTIKAILV STLADI
Source	E.coli
Target Names	arcB
Protein Names	Recommended name: Ornithine carbamoyltransferase, catabolic Short name= OTCase EC= 2.1.3.3
Expression Region	2-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.