



Recombinant Zea mays Zein-alpha GZ19AB11

Product Code	CSB-MP357383ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08416
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	TIFTQCSQA PIASLLPPYL SSAVSSVCEN PILQPYRIQQ AIAAGILPLS PLFLQQSSAL LQQLPLVHLL AQNIRAQQLQ QLVLANLAAY SQQQQFLPFN QLGSLNSASY LQQQQLPFSQ LPAAYPQQFL PFNQLAALNS PAYLQQQQLL PFSQLAGVSP ATFLTQPQLL PFYQHVAPNA GTLLQLQQLL PFNQLALTNP AVFYQQPIIG GALF
Source	Mammalian cell
Protein Names	Recommended name: Zein-alpha GZ19AB11 Alternative name(s): 19 kDa zein GZ19AB11
Expression Region	22-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.