



Recombinant Succinyl-CoA ligase [ADP-forming] subunit alpha (sucD)

Product Code	CSB-MP357582TNR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09143
Product Type	Recombinant Protein
Immunogen Species	Thermus thermophilus
Purity	>85% (SDS-PAGE)
Sequence	MILVNRETRV LVQGITGREG QFHTKQMLDY GTKIVAGVTP GKGGTEVLGV PVYDTVKEAV AHHEVDASII FVPAPAAADA ALEAAHAGIP LIVLITEGIP TLDMVRAVEE IKALGSRLIG GNCPGIISAE ETKIGIMPGH VFKRGRVGII SRSGTLTYEA AAALSQAGLG TTTTVGIGGD PVIGTTFKDL LPLFNEDPET EAVVLIGEIG GSDEEEAAAW VKDHMKKPVV GFIGGRSAPK GKRMGHAGAI IMGNVGTPEK KLRAFAEAGI PVADTIDEIV ELVKKALG
Source	Mammalian cell
Target Names	sucD
Protein Names	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit alpha EC=6.2.1.5 Alternative name(s): Succinyl-CoA synthetase subunit alpha Short name= SCS-alpha
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.