



Recombinant Oreochromis mossambicus Prolactin-1 (prl1)

Product Code	CSB-EP357641OES-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09319
Product Type	Recombinant Protein
Immunogen Species	Oreochromis mossambicus (Mozambique tilapia) (Tilapia mossambica)
Purity	>85% (SDS-PAGE)
Sequence	VPINEL FERASQHSDK LHSLSTTLTQ ELDSHFPPIG RVIMPRPAMC HTSSLQTPID KDQALQVSES DLMSLARSLL QAWS DPLVVL SSSASTLPH AQSTIFNKIQ EMQQYSKSLK DGLDVLSSKM GSPAQAITS L PYRGGTNLGH DKITKLINFN FLLSCLRRDS HKIDSFLKVL RCRAAKMQPE MC
Source	E.coli
Target Names	prl1
Protein Names	Recommended name: Prolactin-1 Alternative name(s): PRL-188 Prolactin I
Expression Region	25-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.