



Recombinant Enterobacteria phage T4 Tail lysozyme (25)

Product Code	CSB-EP357674EDZ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09425
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Purity	>85% (SDS-PAGE)
Sequence	MANINKLYSD IDPEMKMDWN KDVSRSGLR SIKNSLLGII TTRKGSRPFD PEFGCDLSDQ LFENMTPLTA DTVERNIESA VRNYEPRIDK LAVNVIPVYD DYTLIVEIRF SVIDNPDDIE QIKLQLASSN RV
Source	E.coli
Target Names	25
Protein Names	Recommended name: Tail lysozyme EC= 3.2.1.17 Alternative name(s): Outer wedge of baseplate protein Protein Gp25
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.