



Recombinant Phospholipase C (plc)

Product Code	CSB-MP357713BQJ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09598
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus
Purity	≥85% (SDS-PAGE)
Sequence	WS AEDKHKEGVN SHLWIVNRAI DIMSRNTTLV KQDRVAQLNE WRTELENGIY AADYENPYD NSTFASHFYD PDNGKTYIPF AKQAKETGAK YFKLAGESYK NKDMKQAFFY LGLSLHYLGD VNQPMHAANF TNLSYPQGFH SKYENFVDTI KDNYKVT DGN GYWNWKG TNP EEWIHGAAVV AKQDYS GIVN DNTKDW FVKA AVSQEYADKW RAEVTPMTGK RLMDAQRVTA GYIQLWFD TY GDR
Source	Mammalian cell
Target Names	plc
Protein Names	Recommended name: Phospholipase C Short name= PLC EC= 3.1.4.3 Alternative name(s): Cereolysin A Phosphatidylcholine cholinephosphohydrolase
Expression Region	39-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.