



Recombinant *Saccharomyces cerevisiae* Tubulin alpha-3 chain (TUB3)

Product Code	CSB-EP357747SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09734
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MREVISINVG QAGCQIGNAC WELYSLEHGI KEDGHLEDGL SKPKGGEEGF STFFHETGYG KFPRAIYVD LEPNVIDEVR TGRFKELFHP EQLINGKEDA ANNYARGHYT VGREIVDEVE ERIKMDADC DGLQGFLFTH SLGGGTGSGL GSLLENLSY EYGKSKLEF AVYPAPQLST SVVEPYNTVL TTHTTLEHAD CTFMVDNEAI YDICKRNLGI SRPSFSNLNG LIAQVISSVT ASLRFDGSLN VDLNEFQTNL VPYPRIHFPL VSYAPILSKK RATHESNSVS EITNACFEPG NQMVKCDPTK GKYMANCLLY RGDVVTRDVQ RAVEQVKNKK TVQMVDWCPT GFKIGICYEP PSVIPSELA NVDRVCMLS NNTAIADAWK RIDQKFDLMY AKRAFVHWYV GEGMEEGEFT EAREDLAAL RDIIEVGADS YAEEF
Source	E.coli
Target Names	TUB3
Protein Names	Recommended name: Tubulin alpha-3 chain
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.