



Recombinant Shufflon protein C

Product Code	CSB-MP357752ENL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09749
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	≥85% (SDS-PAGE)
Sequence	MKKYDRGWAS LETGAALLIV MLLIAWGAGI WQDYIQTKGW QTEARLVSNW TSAARSYIGK NYTTLQGSST TTPAVITTT MLKNTGFLSS GFTETNSEGQ RLQAYVVRNA QNPELLQAMV VSSGGTPYPV KALIQMAKDI TTGLGGYIQD GKTATGALRS WSVALSNYGA KSGNGHIAVL LSTDELSGAA EDTDRLYRFQ VNGRPDLNKM HTAIDMGSSN LNNVGAVNAQ TGNFSGNVNG VNGTFSGQVK GNSGNFDVNV TAGGDIRSNN GWLITRNSKG WLNETHGGGF YMSDGSWVRS VNNKGIYTGG QVKGGTVRAD GRLYTGEYLQ LERTAVAGAS CSPNGLVGRD NTGAILSCQS GTWGAPKIQF TTQTYNLAKN TRNLRLGVHA YCSWTYLNLS PFGGFQQVYS DQNNVWYVSN YAWGNYESGG TISVTCLNLP GAGA
Source	Mammalian cell
Protein Names	Recommended name: Shufflon protein C
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.